URGENT FIELD SAFETY NOTICE

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Sure-Pro® and Sure-Pro Ultra® Embryo Replacement Catheters

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Smiths Medical is providing this Urgent Field Safety Notice to advise its customers of a Field Safety Corrective Action for Sure-Pro® and Sure-Pro Ultra® Embryo Replacement Catheters ("Sure-Pro Catheters"). Smiths Medical is voluntarily taking this Action with the knowledge of the relevant Regulatory Agencies.

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Smiths Medical has become aware of a small number of users encountering resistance when passing the inner catheter through the outer sheath. A limited number of Sure-Pro Catheters were manufactured with the tip of the inner catheter slightly larger than the inner diameter of the outer sheath. This issue does not involve every device listed on the list of affected products; however all products listed have the potential to exhibit the issue.

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Smiths Medical has received no reports of patient injury related to this issue.

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Advice on Action to be Taken by the User:

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Subject to this Urgent Field Safety Notice, Smiths Medical is requiring its customers to inspect their inventory and choose one of the options below:

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1. Remove the inner catheter and outer sheath from the pre-packed blister tray as shown below.

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
2. Gently insert the inner catheter into the outer sheath and progress the inner catheter along the full length of the outer sheath as shown below

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3. Continue to progress the inner catheter completely through the length of the outer sheath until the tip of inner catheter is visibly protruding through the end of the outer sheath as shown below. If the inner catheter slides through the outer sheath with no significant resistance, proceed with use of the product pursuant to the Instructions For Use. If significant resistance is felt as the inner catheter protrudes through the outer sheath, discard the device and continue the procedure with an alternative product.

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4. Complete and return the Confirmation Form (Attachment 2), by fax at +44 (0)1233 722153 or by email to SurePro@smiths-medical.com, to acknowledge receipt and understanding of the Field Safety Notice and to provide details of affected product in inventory.

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Option 2 – Return Unused Affected Product for Replacement or Credit

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1. If you choose not to perform the inspection as described above, remove all affected products from use.

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2. Complete and return the Confirmation Form (Attachment 2), by fax at +44 (0)1233 722153 or by email to SurePro@smiths-medical.com, to acknowledge receipt and understanding of the Field Safety Notice and to provide details of affected product in inventory. Smiths Medical will contact you to arrange the return of unused affected product and arrange credit or replacement as requested.

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Advice on Action to be Taken by the Distributor:

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1. Inspect your inventory and quarantine all affected products.

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2. Perform a count of affected products currently in inventory and complete and return the attached Confirmation Form (Attachment 2) by fax to +44 (0) 1233 722153 or by email to SurePro@smiths-medical.com. Following receipt of this form, Smiths Medical will contact you to arrange return of the identified affected products.

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3. Send a copy of this Urgent Field Safety Notice to your customers identified as having received affected product.

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4. Provide email confirmation of your customer notifications to Smiths Medical at SurePro@smiths-medical.com.

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Transmission of this Urgent Field Safety Notice

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This notice needs to be passed on to all personnel who need to be aware within your organization, including points of use or to any organization where the potentially affected devices have been transferred. If you or your facility has distributed these affected products to other persons or facilities, please promptly forward the recipients a copy of this Urgent Field Safety Notice.

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Please maintain awareness of this Notice and resulting action for an appropriate period to ensure effectiveness of this Correction

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Customers should report any issues with these products to Smiths Medical's Global Complaint Department at +00 800 76 48 47 00 or globalcomplaints@smiths-medical.com.

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If you should have any questions regarding this information, please contact Smiths Medical's Customer Service Department at +44 (0)1233 722267

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Smiths Medical is committed to providing quality products and service to its customers. We apologise for any inconvenience this situation may have caused.

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Sincerely,

Michael Herbert
Regional Director, Quality Systems
Smiths Medical International

Enclosures: Attachment 1 – Affected Product Reorder Codes and Lot Numbers
 Attachment 2 – Urgent Field Safety Notice Confirmation Form

ATTACHMENT 1:

Příloha č. 1:

Affected Product Reorder Codes and Lot Numbers

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Smiths Medical International
Boundary Road,
Hythe,
Kent, CT21 6JL
United Kingdom

URGENT FIELD SAFETY NOTICE CONFIRMATION FORM

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Customer Identification No. _____

IČ zákazníka

Please complete and return this Form by fax to +44 (0)1233 722153 or by sending an electronic copy via email to SurePro@smiths-medical.com

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