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	<b>SUPPLEMENTAL NOTES FOR USE OF STERILE BARE FIBERS</b>	Code: 96039
		Version: 1 / 1-2273
		Valid from: 24.02.2015

**A. Follow the instructions**

Before attempting to use a single use or a reusable sterile bare fiber please read carefully the accompanying Operating Manual (Instructions for use, User Manual) for that fiber.

**B. Checking the fiber**

Please note that it is of critical importance to check before every use the optical fiber for laser light leakage and shape of the target (pilot) beam as described in the Operating Manual (Instructions for use, User Manual) for that fiber. Any discrepancy may indicate that the fiber is mechanically damaged and not suitable for use.

 **WARNING!**  
When a damaged fiber is used this may result in overheating of the fiber and may cause harm to the patient.

**C. Preparing (re-processing) of re-usable fibers**

Please note that it is of critical importance to strictly follow the procedure for preparing (re-processing) of re-usable fibers as described in the Operating Manual (Instructions for use, User Manual) for that fiber.

 **WARNING!**  
Failure to follow the fiber preparation (re-processing) procedure may result in mechanical damage to the fiber.

For example, when the Operating Manual (Instructions for use, User Manual) instructs the user to strip the fiber jacket by a certain length, it should be performed in a single swift motion. Do not attempt to strip the fiber in multiple pulls. In case the stripping was not successful, cut off the unsuccessfully stripped end of the fiber and repeat the complete preparation procedure.

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